BLOOD PRESSURE MONITOR | Talking

Instruction Manual
Model: BD959

Go to biosmedical.com
for a **FREE** Log Book!

BIOS VP of Marketing - Mark Beaton accepts the 2017 Hypertension Canada Certificate of Excellence from Angelique Berg, CEO of Hypertension Canada

Trusted by Canadians for 3 Generations

At BIOS Diagnostics™, we are proud of our legacy in blood pressure monitoring in Canada. From the early 1930's to 1987 we manufactured "Tycos" brand professional blood pressure equipment for doctors and hospitals in Canada.

In the 1970's we pioneered the first blood pressure devices for monitoring at home, and in the 1980's we introduced digital technology in Canada. We haven't been counting, but we know that millions of our home-use monitors have been used by Canadians in the last 30 years.

All BIOS Diagnostics™ devices are developed in collaboration with physicians and clinically tested to prove their measurement accuracy. For more information on clinical tests and other BIOS medical products, visit our website at **www.biosmedical.com**.

If you have questions about this device or blood pressure monitoring at home, email us at **support@biosmedical.com** or call the **BIOS Medical Hotline 1-866-536-2289**.

Blood Pressure Monitor – Talking Instruction Manual

Table of Contents

1. Introduction

- 1.1 Features
- 1.2 Important Information
 - 1.2A Safety Information
 - 1.2B Care of the Device
 - 1.2C Comparing Readings to Other Blood Pressure Devices
 - 1.2D Calibration
- 1.3 What Do Your Numbers Mean?
- 1.4 Normal Blood Pressure Values
- 1.5 Common Blood Pressure Questions and Answers

2. Getting Started

- 2.1 About the BD959
- 2.2 About the LCD Screen
- 2.3 Powering the Device

3. Using the Device

- 3.1 Setting Date, and Time
- 3.2 Setting the Talking Language
- 3.3 Setting Voice Volume
- 3.4 Obtaining Accurate Measurements
 - 3.4A Tips on Taking Accurate Measurements
 - 3.4B Common Sources of Error
 - 3.4C Fitting the Wide Range Cuff
- 3.4 Measuring Your Blood Pressure
- 3.5 Discontinuing a Measurement
- 3.6 Determining Your “Real” Average Blood Pressure At Home
- 3.7 Recalling the Memory
- 3.8 Clear Measurements from Memory
- 3.9 Blood Pressure Assessment Indicator

4. Error Messages/Malfunctions

5. Care and Maintenance

6. 4 Year Limited Warranty

7. Technical Specifications

1. Introduction

Thank you for purchasing the BIOS Diagnostics™ Blood Pressure Monitor - Talking. Designed for convenient and easy operation, this device provides a simple, yet accurate method to measure your blood pressure.

Your blood pressure is an important parameter that can be used to monitor your health. This device enables you to monitor your blood pressure regularly, and maintain a record of your blood pressure measurements. You can then use this record to assist your physician in diagnosing and maintaining a healthy blood pressure level.

Indications for Use:

This blood pressure monitor (BD959) is intended to measure the systolic/diastolic pressure, and pulse rate of an adult individual (over the age of 12) by using a non-invasive oscillometric technology. This device is designed to be portable, and used in both home and professional environments for every day blood pressure monitoring. This device is not suitable for neonate, pregnancy, or pre-eclampsia.

1.1 Features

The BD959 uses oscillometric technology to measure the arterial blood pressure and pulse rate. The cuff is wrapped around the upper arm, and automatically inflated by the air pump. The device senses weak fluctuations of the pressure in the cuff produced by extension and contraction of the artery within the arm in response to each heartbeat. The amplitude of the pressure waves is measured, converted in millimeters of the mercury column, and is shown on the display.

- **Memory Features** 1 user, 199 blood pressure readings with time and date.
- **Audio Guided Measuring & Memory** choose from English or French.
- **BP Assessment Indicator** displays the range between which your blood pressure values lie, according to Hypertension Canada & the World Health Organization (WHO).
- **Detects Irregular Heartbeat**
- **BIOS Averaging** - Averages the last 3 readings.

This device is easy to use and has been proven in clinical studies to provide excellent accuracy. Before using the BD959, read this instruction manual carefully and keep it in a safe place.

The blood pressure monitor meets the accuracy requirements of Hypertension Canada and has been tested for clinical accuracy.

1.2 Important Information

Refer to the following sections to learn about important safety instructions and how to take care of the BIOS Diagnostics™ Blood Pressure Monitor - Talking.

1.2A Safety Information

- Self-measurement means control, not diagnosis or treatment. Your values must always be discussed with your doctor or a physician who is familiar with your family history.
- If you are undergoing medical treatment and receiving medication, consult your doctor to determine the most appropriate time to measure your blood pressure. Never alter the dosages of any medication without direction from your doctor.
- Your blood pressure depends on several factors, such as age, gender, weight and physical condition. It also depends on the environment and your state of mind at the time of measurement. In general, your blood pressure is lower when you are asleep and higher when you are active. Your blood pressure may be higher when recorded at a hospital or a clinic and may be lower when measured in the relaxing comfort of your home. Due to these variations, we recommend that you record your blood pressure regularly at home as well as at your doctor's clinic.
- Try to record your blood pressure regularly at the same time of the day and under the same conditions. This will help

your physician detect any extreme variations in your blood pressure and thus treat you accordingly.

- Morning Hypertension (>130/80 mmHg): Recently, several studies have identified elevated cardiovascular risks (heart failure, stroke, angina) associated with "morning hypertension". There is a typical rise in blood pressure during the physiological changes from sleep to arising for the day.
- The ideal time to measure your blood pressure is in the morning just after you wake up, before breakfast and any physical activity, and in the absence of the urge to urinate. If this is not possible, try to take the measurements later in the morning, before you start any physical activity. Relax for a few minutes before you record your blood pressure.
- The cuff should not be applied over wounds as it may cause further injury.
- **Do not** attach the cuff to an arm being used for IV infusions or any other intravascular access, therapy or an arteriovenous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the individual.
- Your blood pressure increases or decreases under the following circumstances:
 - Blood pressure is higher than normal:
 - When you are excited, nervous, or tense
 - While taking a bath
 - During and after exercise or strenuous physical activity
 - When it is cold
 - Within one hour after meals
 - After drinking tea, coffee, or other caffeinated drinks
 - After smoking tobacco
 - When your bladder is full
 - Blood pressure is lower than normal:
 - After consuming alcohol
 - After taking a bath
- The pulse display is not suitable for checking the frequency of heart pacemakers.
- If you have been diagnosed with a severe arrhythmia or irregular heartbeat, vascular constriction, liver disorders, or diabetes, have a cardiac pacemaker, or are pregnant, measurements made with this instrument should only be evaluated after consultation with your doctor.
- Take care while handling the batteries in the device. Incorrect usage may cause battery fluid leakage. To prevent such accidents, refer to the following instructions:
 - Insert batteries with the correct polarity.
 - Turn off power after use. Remove and store the batteries if you are not planning to use the device for an extended period of time.
 - **Do not** mix different types, brands, or size of batteries. This may cause damage to the product
 - **Do not** mix old and new batteries.
 - Remove batteries and dispose of them according to the proper regulations in your area.
 - **Do not** disassemble batteries or expose them to heat or fire.
 - **Do not** short-circuit the batteries.
 - **Do not** use rechargeable batteries.

1.2B Care of the Device

For prolonged life of your blood pressure monitor, note the following instructions:

- **Do not** drop or bang the unit. Prevent sudden jerks, jars, or shocks to the device to prevent damage.
- **Do not** insert any foreign objects in any device openings or vents.
- **Do not** disassemble the unit.
- If the unit has been stored at very low or freezing temperatures, allow to reach room temperature before using it.

- **Do not** store the unit in direct sunlight, high humidity, or in places with a lot of dust.
- Clean the device with a soft dry cloth. **Do not** use gasoline, thinner or similar solvents. Carefully remove spots on the cuff with a damp cloth and soap. **Do not** wash the cuff.
- **Do not** use the device if you think it is damaged or if anything appears unusual.
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed.
- Using the unit in the immediate vicinity of mobile phones, microwave appliances, or other devices with strong electromagnetic fields may result in impaired functioning.
- **Do not** use this device close to strong electromagnetic fields, such as mobile telephones or radio installations.
- Keep a distance from such devices when using this unit.

1.2C Comparing Readings to Other Blood Pressure Devices

Many questions arise when two blood pressure devices are compared in an effort to check accuracy. An accurate comparison requires repeatable measurements under the same conditions to a “reference device” with known accuracy. Significant time is required to reduce naturally occurring blood pressure variability during the test. The subject should be seated comfortably with feet flat on the floor, and have rested for 5 minutes before the first reading to allow blood pressure levels to stabilize. The patient's back, elbow, and forearm should be supported, and the middle of the cuff should be at the level of the right atrium. There should be no talking or moving during the measurement and if comparing to an aneroid gauge or mercury column, observers should avoid parallax and be careful not to round measurements.

The most accurate way to compare devices is to take two readings at the same time. However, most people and doctor's offices do not have the equipment necessary to measure blood pressure from two devices simultaneously. To take sequential measurements properly requires a pair of initial measurements to determine the subject's blood pressure level: first with the reference equipment, followed by 60 seconds, then with the monitor-under-test. The actual accuracy test requires three pairs of measurements with 60 seconds between measurements. These measurements are averaged and a comparison can be made. Since most people tend to relax and their blood pressure falls with subsequent measurements, following this protocol reduces these natural changes in blood pressure levels. The standard technical error of both consumer and professional devices is normally ± 3 mmHg, so a discrepancy of 6 mmHg is acceptable even when the devices are working within their specifications.

Any comparisons without a known “reference device” and not following the procedures described above will yield unreliable results. In addition, to do an accuracy test properly the reference device must also be tested to a known reference to confirm its accuracy, prior to being used as the reference for comparisons.

1.2D Calibration

Digital blood pressure monitors do not require regular recalibration, unless the product has been dropped and internal parts have been damaged. If the unit turns on and does not display an error code, the product is working properly. In extremely rare cases, the cuff may have developed a pin-hole leak, or the gasket where the cuff connector enters the monitor may not have a proper seal; both of these leaking air issues will potentially cause errors in accuracy, but otherwise the product will work accurately without drifting out of calibration.

1.3 What do your Numbers Mean?

Blood pressure is the pressure in your blood vessels while blood circulates throughout your body. High blood pressure or “Hypertension” is the pressure at which one's normal average blood pressure is considered too high and other health risks including: heart attack, stroke, dementia, kidney failure, heart disease and erectile dysfunction may occur. It is expressed as two numbers: systolic/diastolic 120 mmHg/80 mmHg (mmHg= millimeters of mercury). “Systolic” numbers refer to the pressure on the walls of your arteries while the heart is contracting and pushing blood. “Diastolic” pressure is the lower number when the heart is at rest and relaxed. A simple way to understand this is to picture a garden hose. When the tap is turned on, the immediate pressure on the walls of the hose is the “systolic” number, and when the tap is turned off it is the “diastolic” number.

There are many different causes of high blood pressure. We differentiate between common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.

1.4 Normal Blood Pressure Values

When rested, and measuring your blood pressure at home, measurements with diastolic pressure readings above 80 mmHg or systolic readings over 130 mmHg are considered high. If you obtain consistent readings in this range, consult your doctor immediately. High blood pressure values over time can damage blood vessels, vital organs such as the kidney, and your heart. With blood pressure values that are too low (i.e., systolic values under 105 mmHg or diastolic values under 60 mmHg), consult with your doctor.

Systolic	Diastolic	Comment
Below 120	Less than 80	This range is considered "Normal" and ideal.
120 - 129	70 - 79	This range is considered "Elevated" or "Pre-hypertension" . Discuss with your health care professional. Lifestyle modifications may be required to avoid advancing into hypertension.
130 - 139	80 - 89	This is in the "Hypertension" range. Discuss with your health care professional. Medication(s) and lifestyle modifications are typical treatments.
140 - 179	90 -119	Discuss with your medical professional, medication(s) and lifestyle modifications are necessary to control your hypertension.
180 and above	120 and above	This is in the "Hypertensive Crisis" range. Consult your doctor immediately, or go to the hospital if you are experiencing any chest pain, shortness of breath, or symptoms of stroke.

Adapted to include new information as outlined in Hypertension Canada's 2025 Understanding and measuring your blood pressure guideline.

Note: A diagnosis of high blood pressure must be confirmed with a medical professional. A doctor should evaluate any unusual blood pressure readings. Additionally, lower targets may be appropriate for some populations such as African-Americans, the elderly, or patients with underlying issues such as diabetes mellitus or chronic kidney disease.

Important for Canadians:


*** Hypertension measured at home $\geq$ 130/80**

For further information, see our website www.biosmedical.com.

1.5 Common Blood Pressure Questions and Answers

a) Why is my blood pressure reading always different?

Your blood pressure changes constantly. It is quite normal for blood pressure to fluctuate as much as 50 mmHg throughout the day. Blood pressure is normally lowest at night, but increases during waking hours when the stress and activities of everyday life are highest.

b) Why is the doctor's reading different from the reading taken at home?

Your blood pressure can vary due to the environment (temperature, nervous condition). When measuring blood pressure at the doctor's office, it is possible for blood pressure to increase due to anxiety and tension, this is known as "White Coat Hypertension". Home blood pressure monitoring has been known to predict health outcomes better than office blood pressure measurements.


c) Why should I monitor blood pressure at home?

One or two readings will not provide a true indication of your normal blood pressure. It is important to take regular, daily measurements and to keep records over a period of time. This information can be used to assist your physician in diagnosing and preventing potential health problems. Expert opinion suggests patients with normal blood pressure can benefit from one week of measuring every three months.

2. Getting Started

2.1 About the BD959

This section describes the various components of the Blood Pressure Monitor – Talking.


Upper arm cuff:


Wide range cuff for arm circumference 22-45 cm (8.7"-16.5")

2.2 About the LCD Screen

The LCD screen displays the systolic and diastolic blood pressure measurements along with your heart rate. It also displays previously recorded measurements and the date and time, when the appropriate button is pressed.

- | | |
|--------------------------------|--------------------------------|
| 1. Time & Date Display | 7. Irregular Heartbeat Icon |
| 2. Voice Icon | 8. Memory Icon |
| 3. Battery Power Icon | 9. Pulse Display/Memory Number |
| 4. Inflation & Deflation Icons | 10. Diastolic Pressure (DIA) |
| 5. Cuff Wrap Icon | 11. Systolic Pressure (SYS) |
| 6. Movement Error Icon | 12. BP Assessment Indicator |


	If the cuff is wrapped incorrectly, it may cause unreliable values. When the  is displayed, the cuff is tight enough. When the  is displayed, the cuff is too loose, and needs to be re-applied.
	Appears when your body moves during a measurement. Remove the arm cuff, wait 2-3 minutes, sit still, and try again.
	Appears when an irregular rhythm is detected during a measurement. If it continues to appear, consulting with your physician is recommended.
	Heartbeat Icon flashes while taking a measurement.

2.3 Powering the Device

If using the USB-C power cable and adapter (included), follow these steps to power the device on.


1. Connect the USB-C cable head to the USB-C cable port on the blood pressure device.
2. Connect the opposite cable end to the wall adapter.
3. Plug the wall adapter into the socket.

Note: The device does not hold a charge. It must remain connected to the power socket to function.

If using batteries (not included), follow these steps to insert four "AA" batteries into the device.

1. Open the battery compartment cover in the direction shown.
2. Insert four "AA" batteries ensuring the correct polarity is followed.
3. Replace the battery compartment cover.

3. Using the Device

3.1 Setting Date, and Time

The monitor stores measuring results automatically after a reading is taken. For accurate memory information, it is necessary to reset the time and date whenever a new power source is introduced (i.e. new batteries are inserted, or power supply through adapter was disconnected/reconnected).

1. With the power off, press and hold both the **Start/Stop button** , and the **Memory button**  for approximately 3 seconds until the year begins to flash.
2. Use the **Memory button**  to select the year, and the **Start/Stop button**  to save selection.
3. Use the same button selection method to set the Month, and Date, followed by the Time (Hour, then minute). (**Memory button**  to select, **Start/Stop button**  to save selection).

3.2 Setting the Talking Language

After setting the Time, the next selection will be Language.

1. Toggle between En (English) and Fr (French) using the **Memory button** . Use the **Start/Stop button**  to save selection.

3.3 Setting Voice Volume

Once the Language setting has been saved, the final setting is Voice Volume.

- 1 Using the **Memory button** , cycle through the voice volume levels to hear an example of the sound output.
2. Then select and save the choice using the **Start/Stop button** .

There are 5 volume selection levels: 1, 2, 3, 4, and Off

Note: It is by updating this selection setting that you can turn the voice off to make the device silent, not within the Talking Language setting.

3.4 Obtaining Accurate Measurements

Your blood pressure can vary based on numerous factors, physiological conditions, and your surroundings. Follow these guidelines to obtain accurate and error-free measurements of your blood pressure and pulse rate.

3.4A Tips on Taking Accurate Measurements

In the morning before breakfast, 2 hours after dinner, before taking medication.

Empty bladder (if necessary).

Avoid coffee and smoking within the hour, and no exercise 30 minutes before measuring.

Rest quietly for 5 minutes. Remain calm and quiet while the measurement is in process.

Do not speak while taking the measurement.

Take measurements on the non dominant arm.

Sit with legs uncrossed so as not to restrict blood flow.

Sit with back supported and measurement arm resting on a table. Sit with feet flat on the floor.

Ensure that the cuff is level with the heart while the arm is supported on the table.

3.4B Common Sources of Error

All efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.

ATTENTION!


Comparable blood pressure measurements always require the same conditions with a peaceful and calm environment. Ensure that you take measurements under the same conditions to obtain accurate and reliable readings.

- If the arm artery lies considerably lower or higher than the heart, an erroneous value of blood pressure is measured. Each 15 cm difference in height results in a measurement error of 10 mmHg.
- A loose cuff causes false measurement values.
- With repeated measurements, blood accumulates in the arm, which can lead to false results. Consecutive blood pressure measurements should be repeated after at least a 45 second pause or after the arm has been held up in order to allow the accumulated blood to flow away.


3.4C Fitting the Wide Range Cuff

NOTE: Visit www.biosmedical.com for a video that illustrates how to properly fit the cuff.


- Putting the cuff Velcro® side down on a table, the words on the cuff should be facing right side up. Put the end of the cuff (with fastener) through the metal ring, making the cuff a cylinder. (Ignore this step if your cuff is already set up.) Proper assembly allows the Velcro® to match up properly. Measure on your non-dominant arm, unless there is a >10 mmHg difference with the other arm, in which case use the arm with the higher pressure.


- b) Place the cuff around your bare arm. Make sure the bottom edge of the cuff is about 1" (2-3 cm) above the elbow joint. Adjust the cuff so that the rubber tubing under the cuff lies over the brachial artery, which runs on the inside of the arm (see Fig. B). The red material should be over the brachial artery.


- c) Pull the cuff and tighten it by attaching the Velcro® fastener.


- d) The cuff should fit snugly around the bare arm, but not too tight. You should be able to fit two fingers under the cuff. If the cuff is the wrong size, the device will not measure your blood pressure accurately. Contact the store, or BIOS Medical for other sized cuffs.

- e) Place the arm on the table (palm facing upwards) so that the cuff is at the same level as the heart. Make sure there is no kink in the hose.

- f) You can adjust the level of your arm by putting a cushion under your arm. The cuff should be at heart level.

- g) Remain seated in a comfortable room temperature for at least 5 minutes, f) then start the measurement.

- h) For those who cannot put the cuff on the left arm, put it on the right arm as shown.

- i) More than 6 consecutive measurements will cause blood accumulation in the lower arm which will affect the measuring results. To improve reading accuracy, raise the arm being measured, squeeze and relax your hand several times, then take another measurement. Another option is to take the cuff off and wait at least 5 minutes before repeating measurement.

i)

- j) If this device was stored in low temperature, it is necessary to leave it in room temperature for at least 1 hour, otherwise the measurement can be inaccurate.

Comment:

Continue to use the same arm for comparisons. It is not unusual for there to be a difference in blood pressure between arms. Initially, check BP on both arms. If one arm is 10mmHg higher, continue to measure on that arm.


Comparable blood pressure measurements always require the same conditions (relax for several minutes before taking a measurement).

ATTENTION: Do not use a cuff other than the original cuff contained in this kit, or one purchased from the manufacturer!


3.5 Measuring Your Blood Pressure

After the cuff has been appropriately positioned, the measurement can begin:

- Remain seated in a comfortable room temperature location for at least 5 minutes, and then start the measurement.
- Press the **Start/Stop button**  to start the measurement. The device will quickly complete a self-diagnosis check where all icons on the screen will appear, followed by the deflation icon  flickering, ensuring the device is prepared for the reading. Once completed the inflation icon  will flicker, and the device will begin to inflate the cuff. Measurement procedure has begun.
- The monitor starts measuring automatically after the inflation has finished. The heartbeat icon  flickers, and the pressure release value decreases gradually.


Full screen display state

The zeroing state


The inflation state

The measuring state

- When the measurement has been completed, the systolic, diastolic, and pulse rate will appear on the display.

- The measurement reading remains on the display until you switch off the device. If no button is pressed, the device switches itself off in approximately 3 minutes.

3.6 Discontinuing a Measurement

You can stop the measurement at any time by pressing the **Start/Stop button**  (e.g. if you feel uneasy or an unpleasant pressure sensation). The device will immediately decrease the cuff pressure automatically.

3.7 Determining Your “Real” Average Blood Pressure at Home

It is normal for blood pressure to vary significantly in the middle of the day when most people are busy with their daily tasks. Hypertension Canada recommends measuring in the morning and evening to avoid variability. At the end of a measurement, this instrument automatically stores each result, including date and time.

AM

Take: 2 measurements, 1 minute apart

- Empty bladder (If necessary)
- In the morning before breakfast, before taking medication
- Sit with back supported and measurement arm resting on a table. Sit with feet flat on the floor.

PM

Take: 2 measurements, 1 minute apart

- Empty bladder (If necessary)
- 2 hours after dinner, before taking medication
- Sit with back supported and measurement arm resting on a table. Sit with feet flat on the floor.
- Avoid coffee and smoking within the hour, and no exercise 30 minutes before measuring.

RESULTS

Discard Day 1 measurements.
Average your day 2-7 measurements

= Average

< 130/80 mmHg = NO

No Hypertension

≥ 130/80 mmHg = YES


Yes Hypertension

***Note: If the result is “borderline” repeat the series for confirmation. This data can be used by a medical professional to make a diagnosis of hypertension.**

3.8 Recalling the Memory

This blood pressure monitor automatically stores 199 measurements in its memory. The oldest record will be replaced by the latest measurement value when there are more than 199 sets.


a) Press the **Memory button** (M). The average of the last 3 measurements will be shown.

The average value

b) Continue to press the **Memory button** (M) to cycle through individual measurement history with date and time information (1 through 199).

Note: When reviewing memory, the date and time will show alternatively. Remain on the reading to see both.

The 1st memory


c) To cycle through the memories quickly to reach a specific number, press and hold the **Memory button** (M), releasing when you've reached the desired range.

d) When done reviewing memories, press the **Start/Stop button** (power icon) to shut the monitor off.

3.9 Clear Measurements from Memory

If you are sure that you want to permanently remove all stored memories, enter into the device memory by pressing the **Memory button** (M).

Next, press and hold both the **Start/Stop button** (power icon) and the **Memory button** (M) until the display shows the below information:

Press the **Start/Stop button**  to shut the monitor off.

3.10 Blood Pressure Assessment Indicator

The bars on the left hand side of the display show you the range within which the indicated blood pressure values lies. Depending on the height of the bar, the readout value is either within the normal (green), borderline (yellow) or hypertensive (red) range.

The classification is based on standards adopted from the WHO*
Refer to the charts below for details of the classification.

	SYS (mmHg)	DIA	
Red	180 ▲	120 ▲	3 — Hypertensive Crisis
Red	140-179	90-119	2 — Moderate Hypertension
Red	130-139	80-89	1 — Mild Hypertension
Yellow	120-129	70-79	— Elevated
Green	110-119	70-79	N — Normal
Green	▼110	▼70	L — Optimal

Adapted to include new information as outlined in Hypertension Canada's 2025 Understanding and measuring your blood pressure guideline.

- If the bar lines up with the first or second green bar, your measurement is in the green zone, or "Normal".
- If the bar lines up with the third bar, it is in the yellow zone, and is considered elevated, or pre-hypertension.
- If the bar lines up with the fourth bar, it is in the red "Stage 1 Hypertension" zone.
- If the bar lines up with the fifth bar, it is in the red "Stage 2 Hypertension" zone.
- If the bar lines up with the sixth bar, it is in the red "Hypertensive Crisis" zone.

4. Error Messages / Malfunctions

If an error occurs during a measurement the LCD displays the corresponding error code.



Error	Possible Cause	Remedy
Err 4	Abnormality in pulse wave detection	Ensure that the cuff is fastened properly before starting measurement.
Err 5	Pressurizing error	Check if the cuff connector is correctly installed, or check if there is air leakage from the cuff.
Err 6	Measurement errors caused by arm or body movement	Ensure that both the body and arm are held still for the duration of the measurement.
Err 7	Cuff is loose	Ensure that the cuff is fastened properly before starting measurement.
Err 8	Pressure exceeds the maximum value (300mmHg)	Ensure that the cuff is fastened properly before starting measurement.
	Battery is low	Replace batteries.
	Battery is depleted	Replace batteries.

If problems occur when using the device the following points should be checked, and if necessary, the corresponding measures should be taken.

Malfunction	Remedy
The machine does not work, even after pressing the Start/Stop button, with batteries installed.	1. Check batteries for correct polarity. 2. Check to see if batteries need to be replaced.
The device frequently fails to measure the blood pressure values, or the values measured are too low or too high.	1. Check the positioning of the cuff and cuff tension. The cuff should be snug but not tight. 2. Measure the blood pressure again, ensuring that you have remained motionless for a sufficient amount of time to ensure an accurate reading.
Device pump is working, but the pressure in the cuff is not rising.	Check that the cuff connector is properly installed into the side of the monitor.
Every measurement produces varying results although the instrument functions normally and the values displayed are normal.	Note that blood pressure fluctuates continuously; therefore measurements will show some variability.
Blood pressure values measured differ from those measured by the doctor.	Record daily measurements for consultation with your doctor. Note: Individuals visiting their doctor frequently experience anxiety, which can result in a higher blood pressure reading than at home.

When the determined blood pressure is outside the device's rated range, there is a visual alarm signal on the display screen.


When an alarm signal is triggered, no measurement results will be displayed. Please refer to below table for details.

Alarm Indication	Display Contents and Causes
	① When SYS display area shows "HI", the systolic pressure result exceeds 260 mmHg. ② When DIA display area shows "HI", the diastolic pressure result exceeds 210 mmHg.
	① When SYS display area shows "LO", the systolic pressure result is below 40mmHg. ② When DIA display area shows "LO", the diastolic pressure result is below 20mmHg.
If the above alarm indications show again after resting, and taking a new measurement, consult with your doctor.	

For assistance call BIOS Medical Blood Pressure Hotline: 1-866-536-2289

5. Care and Maintenance

- Do not** expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of stress through twisting or buckling.
- Clean the device with a soft, dry cloth. **Do not** use gas, thinners or similar solvents. Dirt on the cuff should be removed with a damp cloth and soap. **Do not** put in the washing machine or dishwasher. **Do not** submerge in water.
- Handle the tube carefully. **Do not** pull on it. **Do not** allow the tubing to kink and keep it away from sharp edges.
- The monitor contains sensitive parts, and should be treated gently.
- Never open the monitor. This invalidates the manufacturer's warranty.
- Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

6. 4 Year Limited Warranty

This BIOS Diagnostics™ blood pressure monitor has a 4 year limited warranty to be free of defects from the date of purchase. Only applicable to the original owner. This warranty does not include the inflation system including the cuff and inflation bladder. The cuff is warrantied for 2 years. The warranty does not cover damage from misuse or tampering.

100% Satisfaction Guarantee

If at any time, you are not completely satisfied with the performance of this device, call our BIOS Medical Hotline and speak with a customer service person, who will make arrangements to have the device tested or replaced to your full satisfaction.

If you have questions regarding the operation of your monitor call the **BIOS Medical Hotline: 1-866-536-2289**

Should correction be necessary, return the unit with all component pieces. Enclose proof of purchase and \$5.00 for return shipping and insurance.

THERMOR LTD.

Returns Department
16975 Leslie Street
Newmarket, ON L3Y 9A1
www.biosmedical.com
Email: support@biosmedical.com

Please include your name, return address, phone number, and email address. Thermor will test or replace (at Thermor's option) free of charge any parts necessary to correct the defect in material or workmanship.

Please allow 10 days for return shipping.

7. Technical Specifications

Display:	LCD Digital Display
Measurement Method:	Oscillometric Method
Measurement Range:	Diastolic: 20 mmHg-210 mmHg Systolic: 40 mmHg-260 mmHg Cuff pressure: 0 mmHg-300 mmHg
Pulse Rate :	40 beats/min-200 beats/min
Accuracy:	Pressure: within ± 3 mmHg Pulse rate: within 5% of reading value
Operating Environment	Temperature Range: +5°C to +40°C Relative Humidity Range: 15% to 90%, non-condensing Atmospheric Pressure Range: 70kPa to 106kPa.
Storage/Transportation Environment*	-25°C to +5°C, and +5°C to +35°C at a relative humidity up to 90%, non-condensing, >35°C to 70°C at a water vapour pressure up to 5kPa.
Cuff Size	22-45 cm (wide range cuff - included)
Dimensions:	Approx.130 x 96 x 75 (mm)
Weight:	260g (without batteries)
Power Source:	Battery: 4 AA batteries

AC adapter: input 100V-240V~ 50Hz/60Hz
0.35A MAX output 5V  1000mA

Battery Life:

Approximately 300 measurements

Service Life:

Monitor: 5 years (6 measurements a day)
Cuff: 5000 inflations

Reference to Standards:

ISO 13485:2016	ISO 10993-1
ISO 14971:2019	ISO 10993-5
IEC 80601-2-30	ISO 10993-10
IEC 60601-1	ISO 10993-18
IEC 60601-1-2	ISO 10993-23
IEC 60601-1-11	

***Storage/Transportation Environment:**

- 1. When the ambient temperature is 20°C, the time required for the device to warm from the minimum storage temperature (-25°C) until the device is ready for use is 2 hours.*
- 2. When the ambient temperature is 20°C, the time required for the device to cool from the maximum storage temperature (70°C) until the device is ready for use is 2 hours.*

Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use.

Type BF applied part

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

IP20

No protection against ingress of water. Protected against solid foreign objects of 12.5mm and greater.

BIOS I Medical
MANUFACTURED BY/FABRIQUÉ PAR :
THERMOR LTD.
16975 LESLIE STREET
NEWMARKET, ON L3Y 9A1
MADE IN CHINA/FABRIQUÉ EN CHINE
WWW.BIOSMEDICAL.COM
 100% Canadian Owned & Operated
Propriétaire-exploitant canadien à 100 %

Electromagnetic Compatibility Information

This device conforms to IEC 60601-1-2:2014+AMD1:2020 Electromagnetic Compatibility (EMC) standard.

Table 1 – Emission Limits and Compliance

Phenomenon	Emission Limits	Compliance
Conducted and radiated RF Emissions	CISPR 11	Group1, Class B
Harmonic distortion	IEC 61000-3-2	Class A
Voltage fluctuations and flicker	IEC 61000-3-3	Complies

Table 2 – Immunity Test Levels

Phenomenon	Basic EMC Standard	Immunity Test Levels
Electrostatic discharge	IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air for enclosure port
Radiated RF EM fields	IEC 61000-4-3	10 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz for enclosure port
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz and 60 Hz for enclosure port
Proximity magnetic fields	IEC 61000-4-39	See table 4
Electrical fast transients / bursts	IEC 61000-4-4	±2 kV 100 kHz repetition frequency for input a.c. power port
Surges Line-to-line	IEC 61000-4-5	±0.5 kV, ±1 kV for input a.c. power port
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 Vrms 0.15 MHz to 80 MHz 6 Vrms in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz for input a.c. power port
Voltage dips	IEC 61000-4-11	0 % U_T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° for input a.c. power port
		0 % U_T ; 1 cycle and 70 % U_T ; 25/30 cycles single phase: at 0° for input a.c. power port
Voltage interruptions	IEC 61000-4-11	0 % U_T ; 250/300 cycle for input a.c. power port

Table 3 – Test Specifications for Enclosure Port Immunity to RF Wireless Communications device

Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation	Immunity Test Level (V/m)
385	380 to 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	27
450	430 to 470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine	28
710	704 to 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	9
745				
780				
810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	28
870				
930				
1720	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	28
1845				
1970				
2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	28
5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	9
5500				
5785				
If necessary to achieve the immunity test level, the distance between the transmitting antenna and the ME equipment or ME system may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.				
a) For some services, only the uplink frequencies are included.				
b) The carrier shall be modulated using a 50 % duty cycle square wave signal.				
c) As an alternative to FM modulation, the carrier may be pulse modulated using a 50 % duty cycle square wave signal at 18 Hz. While it does not represent actual modulation, it would be worst case.				

Table 4 – Test Specifications for Enclosure Port Immunity to Proximity Magnetic Fields

Test Frequency	Modulation	Immunity Test Level (V/m)
30 kHz	CW	8
134.2 kHz	Pulse modulation ^{a)} 2.1 kHz	65 ^{b)}
13.56 MHz	Pulse modulation ^{a)} 50 kHz	7.5 ^{b)}
a) The carrier shall be modulated using a 50 % duty cycle square wave signal. b) r.m.s., before modulation is applied.		

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