BIOS living™

Massaging Hot/Cold Water Bottle Cover
LF308

INSTRUCTIONS:

When using a battery operated device, basic precautions should always be followed, including the following:

READ ALL INSTRUCTIONS BEFORE USING

- Close supervision is necessary when the unit is used by, on or near children, or persons with disabilities.
- Use this unit only for its intended use as described in this manual.
- **Do not** place or drop into water or any other liquid.
- **Do not** place or store unit where it can fall or be dropped into a tub or sink.
- Turn off unit when not in use.
- Never operate this unit if it has been damaged or if it is not working properly.
- Keep away from heated surfaces.
- This product is a non-professional appliance designed to provide soothing massage to worn muscles. **Do not** use this product as a substitute for medical attention.

SAVE THESE INSTRUCTIONS

Caution: Please read all instructions carefully before operating.

- If you have concerns regarding your health, consult your doctor before using this product.
- Individuals with pacemakers should consult a doctor before use.
- Those with the following conditions should consult their doctor before using this product:
  - Pregnancy
  - Diabetes
  - Heart Disease
  - Malignant Tumor
  - Conditions requiring rest
- If at any time soreness, rash or an increase in pain occurs, discontinue use and consult your physician.
- NEVER use this product directly on swollen or inflamed areas, skin eruptions or an open wound.
- A massage session longer than 15 minutes is not recommended.

BATTERY PRECAUTIONS

- Use only the size and type of batteries specified.
- When installing batteries, observe proper +/- polarities. Incorrect installation of batteries may cause damage to the unit.
- **Do not** mix different types of batteries together (e.g., alkaline with carbon-zinc or old batteries with new ones).
- If the unit is not going to be used for a long period of time, remove batteries to prevent damage due to possible battery leakage.

BATTERY INSTALLATION

1. Unzip the Hot Water Bottle Case at the bottom to reveal the battery compartment.
2. Place 2 "AA" batteries into the battery compartment. (Proper battery orientation is also noted inside the battery compartment housing).

INSTRUCTIONS FOR USE:

TO INSERT WATER BOTTLE:

- Make sure water bottle is empty.
- Pull the velcro apart and put water bottle in through this opening.

TO TAKE OUT THE WATER BOTTLE:

- Make sure water bottle is empty.
- Pull the velcro open and slide water bottle out.

MASSAGE ACTIVATION:

- Massage may be activated by pressing the POWER button. Two massage modes are available.
- Press the POWER button to change modes or turn power ON/OFF.

NOTE: This unit can be used with or without massage at anytime.

MAINTENANCE

TO CLEAN:

- Clean only with a soft, slightly damp cloth. Never allow water or any other liquids to come into contact with the unit.
- **Do not** immerse in any liquid to clean. To clean the fabric cover, wipe with damp cloth, using a mild liquid detergent or suds. Wipe surface with a clean dry cloth.
- Never use abrasive cleaners, brushes, gasoline, kerosene, glass/furniture polish or paint thinner to clean the massager.
- **Do not** attempt to repair this unit. There are no user serviceable parts.

1 YEAR WARRANTY

The BIOS Living™ Massaging Hot/Cold Water Bottle Cover has a 1 year limited warranty to be free of manufacturing defects in materials and workmanship under normal applications for 1 year of the original owner. If this product becomes inoperable due to defect, return the product with all component pieces and proof of purchase to the address listed below. This warranty does not cover any shipping/transport costs. This warranty does not apply if the product is subject to misuse, neglect, rough handling or damage.

Ship the unit prepaid and insured (at owners option) to:

Thermor Ltd.
Returns Department
16975 Leslie Street
Newmarket, ON L3Y 9A1
www.biosmedical.com
Email: support@biosmedical.com
Telephone: 1.866.517.2970

FSC www.fsc.org FPO

Paper Papier

BIOS | Medical
MANUFACTURED BY/FABRIQUÉ PAR :
THERMOR LTD.
16975 LESLIE STREET
NEWMARKET, ON L3Y 9A1
MADE IN CHINA/FABRIQUÉ EN CHINE
WWW.BIOSMEDICAL.COM
100% Canadian Owned & Operated
Propriétaire-exploitant canadien à 100 %